

VUSN GRADUATE COUNCIL LAB COAT ORDER FORM

Orders must be received by July 15, 2016

Name: _____

Address: _____ Phone #: (____) _____

_____ Email: _____

Entering as: _____ Pre-Specialty or ASN-MSN _____ Direct Entry or Post Masters

White Lab Coat w/VUSN embroidery & personalized name - \$55.00 per coat



Women's Staff Length

Five-button closure. Upper left chest pocket and two lower patch pockets. Pass-through pocket openings. One-piece back belt for extra detail. Approximate finished length 36". Imported. 65/35 polyester/cotton, white twill. Home launder. Misses sizes 4-28. *Turn over for sizing information.*



Men's Staff-Length

Four-button closure. Upper left chest pocket and two lower patch pockets. Pass-through pocket openings. Approximate finished length 39". Imported. 65/35 polyester/cotton, white twill. Home launder. Men's regular and tall sizes 34-60. *Turn over for sizing information.*

Women's – 113694-000

Size _____

Quantity _____

Men's – 110528-920

Size _____

Quantity _____

Please print below the name you'd like embroidered on your lab coat (*No certifications are allowed*):
 (*One character per space, max 18 characters*) The font is Roman block, mixed. Please use an initial for your first name if your full name exceeds 18 characters.

--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

Please return your order form to:
VUSN Graduate Council
 c/o Sarah Ramsey
 461 21st Ave. South
 217 Godchaux Hall
 Nashville, TN 37240

of coats ordered _____ x \$55 = \$ _____
 (this includes TN sales tax)

Check enclosed for a total cost = \$ _____

**** Please make checks payable to VUSN ****

Measurements for Sizing a Lab Coat

How to Size Women's Lab Coats

Women's lab coats typically have a different cut than do men's coats. Women buyers should use their dress size to determine their coat size. If you measure yourself, ensure to wrap the tape measure relatively loosely under your armpits, over the shoulder blades, and across the front of your chest.

# Size	4	6	8	10	12	14	16	18	20	22	24	26	28
Bust	34	35	36	37	38	40	41	43	45	48	50	52	54

How to Size Men's Lab Coats

Men should generally use suit jacket measurements to size their lab coats. Look for your size on the inside of a jacket or suit you already own, or you can measure the circumference of your chest. Place the tape under the armpits and wrap it around the shoulder blades in back to account for the widest part of your torso. Height: Students taller than 6-feet may want to opt for a "long" or tall lab coat. The tall designation adds several inches to the length of the sleeves to accommodate longer arms.

Size/Chest	34	36	38	40	42	44	46	48	50	52	54	56	58	60
------------	----	----	----	----	----	----	----	----	----	----	----	----	----	----